






ESwab Collection Instructions:

| Step | Action |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | Perform hand hygiene and put on gloves. |
| 2 | Positively identify the patient using at least two patient identifiers. (i.e. Name, DOB, and/or MRN) |
| 3 | Open the ESwab peel pouch and remove the swab |
| 4 | Collect the patient sample using the swab. Note: Touching the swab applicator below the pink breakpoint should be avoided as it could lead to contamination and incorrect results.  |
| 5 | Remove the screw cap from the ESwab tube and insert the swab all the way to the bottom of the tube.  |
| 6 | While holding the tube away from your face, hold the end of the swab shaft and bend it at a 180-degree angle to break it at the marked breakpoint  |
| 7 | Screw the cap on tightly to prevent leakage.  |
| 8 | Dispose of the swab shaft in a regular trash can |
| 9 | Apply the Specimen information label to the tube. (Do not cover the tube's lot number or expiration date.) Write the date and time of collection on the label.  |
| 10 | Place the tube in a biohazard transport bag |
| 11 | Remove gloves and perform hand hygiene |
| 12 | Promptly transport specimen to the Laboratory. Specimens must be received in laboratory within 24 hours of collection. |

Cautionary Notes per the manufacturer:

- Do not send a dry ESwab as this will lead to unsatisfactory results; Swabs sent with No liquid or very little liquid will be rejected for recollection by the laboratory.
- If the tube spills its contents prior to inserting the swab, the liquid is non-toxic. Put the swab in another tube before sending it to the laboratory and discard the spilled tube.
- If the tube spills after contamination from inserting the swab, follow your facilities instructions for blood and body fluid clean up. Refer to your facility's infection control manual for further instruction.
- If contaminated fluid splashes onto the personnel collecting the sample, treat it as a blood and body fluid exposure. Refer to your facility's infection control manual for further instruction.

- Swabs are the least appropriate specimen for microbiology analysis, as the organisms isolated may only be colonizing the area and are not involved in the infective process.
- Cultures taken from open skin or abscesses may be compromised due to the fact that these lesions are often colonized with a large number of indigenous microbiota.